

Supplementary material

Supplementary Material Appendix 1.

Table 1. All primer and adapter sequences used for the MS-AFLP protocol.

Primer/Adapter Name	Sequence (5' to 3')
EcoRI Adapter Forward	CTCGTATACTGCGTACC
EcoRI Adapter Reverse	AATTGGTACGCAGTA
MspI/HpaII Adaptor Forward	GATCATGAGTCCTGCT
MspI/HpaII Adaptor Reverse	CGAGCAGGACTCATGA
EcoRI Pre-selective Primer	TACTGCGTACCAATTCA
MspI/HpaII Pre-selective Primer	ATCATGAGTCCTGCTCGG
EcoRI Selective Primer	6-FAM-TACTGCGTACCAATTCAGC
EcoRI Selective Primer	5HEX-TACTGCGTACCAATTCACG
MspI Selective Primer	ATCATGAGTCCTGCTCGGTCAT